

Safety and Efficiency of Intravascular Temperature Management for Therapeutic Hypothermia



Donald W. Marion, MD
Consulting Medical Director,
Apothecom Medical Education
301 Edgewater Place, Suite 300
Wakefield, MA

Please address all correspondence to:
Donald W. Marion, MD
35 High Rock Road
Wayland, MA 01778
Email: donmarion1@yahoo.com
Phone: 617 512-0065

Introduction

Therapeutic moderate hypothermia currently is recommended for treatment of patients with out-of-hospital cardiac arrest (OHCA) and for newborns with neonatal hypoxic/ischemic encephalopathy.(2, 23) For patients with OHCA, hypothermia has been shown to result in a 68% greater likelihood of a favorable neurologic outcome early after the event, and a 44% greater likelihood of being alive with a favorable neurologic outcome at 6 months after the arrest, as compared with patients who are not cooled.(1, 7, 26) These studies led the Advanced Life Support Task Force of the International Liaison Committee on Resuscitation (ILCOR) and the American Heart Association (AHA) to recommend that patients with an initial rhythm of ventricular fibrillation following OHCA should be cooled for 12 to 24 hours.(12, 42) For neonates with clinical evidence of hypoxic/ischemic encephalopathy at birth, cooling to 33 -35 C for 48-72 hours, and begun within 6 hours after birth, is recommended based on the results of three separate randomized, controlled trials reported in 2005 by Gluckman et al., Shankaran et al., and Eicher et al..(16, 21, 47)

Therapeutic moderate hypothermia also has been used for the acute care of patients with severe traumatic brain injury in at least nine major clinical trials during the past 15 years.(9, 10, 19, 28, 34, 44, 48, 49, 57) Careful analysis of most of the clinical trials where ICP data is provided suggests that hypothermia is associated with a significant decrease in the incidence of intracranial hypertension.(8, 10, 19, 28, 34, 39, 41, 49, 53-55) And a new metaanalysis completed by the evidence-based guidelines group at the University of Oregon Health Sciences Center, and published in the Third Edition (2007) of the Guidelines for the Management of Severe Traumatic Brain Injury, found that hypothermia was associated with a 46% greater likelihood of a good outcome (GOS of 4-5).(3) Inspired by these findings, a new multicenter prospective randomized clinical trial of therapeutic moderate hypothermia for pediatric traumatic brain injury has been initiated by David Adelson at the University of Pittsburgh. (*R01 NS052478-01A2 (Adelson); Pediatric Traumatic Brain Injury Consortium: Hypothermia*)

But the potential clinical benefits of therapeutic moderate hypothermia for these and other diseases may not be fully realized because of delays in reaching the target temperature.(33) Such delays are inherent with the use of surface cooling techniques, and in most of the adult studies cited above the target temperature was not reached for 6-12 hours after initiation of cooling. Given the unavoidable delays in the discovery, resuscitation, and transport of patients with cardiac arrest (CA), stroke or traumatic brain injury (TBI), it is very likely that significant secondary injury already has occurred by the time hypothermia is initiated. The sooner the patient is cooled to target temperature, the more likely they are to benefit from the therapy. In addition, surface cooling techniques are relatively inefficient, and in their study of cooling using ice packs for patients with OHCA, Holzer reported that they were able to reach target temperature in only 30% of their patients.(1) Study groups for the use of hypothermia in cardiac arrest, neonatal hypoxia/ischemia, and traumatic brain injury all have suggested that reaching the target temperature sooner may allow for hypothermia to be even more beneficial than has already been demonstrated.(2, 3, 23)

Intravascular Cooling Technology

Intravascular temperature management (IVTM) systems have been approved for the rapid cooling of patients treated with therapeutic hypothermia. An early version of IVTM involved cooling catheters inserted into the superior vena cava via the subclavian vein, and a large clinical trial found that this system was 64% more effective than surface cooling techniques for limiting fever in the neuro-intensive care unit.(13) In a subsequent study of patients who were successfully resuscitated after cardiac arrest, Holzer, et al, compared outcomes in 97 patients treated with IVTM for 24 hours to the outcomes in 941 patients who had been treated with "standard" postresuscitation therapy.(27) The likelihood of survival with a favorable neurologic outcome at 30 days after the arrest was 2.15 fold greater in the group receiving IVTM than the standard postresuscitation group.

There now are at least 12 clinical reports of the use of IVTM for the induction and maintenance of therapeutic hypothermia to treat OHCA, acute myocardial infarction (AMI), or ischemic stroke.(Table 1) Among the 8 studies where data is reported, the mean rate of cooling with IVTM is 1.7 C/h.(4, 14, 18, 20, 22, 30, 43,58) Most reports of surface cooling describe a maximum attainable cooling rate of less than 1 C/h. In addition to a more rapid rate of cooling, these studies also show that IVTM allows for tight control of the temperature once the target temperature is reached. In one recent study, IVTM with the Icy catheter (Alsium Corporation) was found to be significantly better than 4 other commonly used temperature management methods for maintaining the target temperature.(58)

And in all but one of the studies, the target temperature could always be reached, but surface cooling methods often were not sufficient to reach the target temperature.(12,17,18,58)

Table 1: Characteristics of endovascular cooling in studies of patients with out-of-hospital cardiac arrest, acute MI, and stroke *

Study	Speed of Cooling	Able to Reach Target Temperature	Maintenance of Target Temperature
Feuchtl (17)	Not Reported	All 19 in IVTM 2/20 SC	Not Reported
Flemming (18)	0.83°C/hr: IVTM 0.09°C/hr: SC	All 31 in IVTM 4/49 in SC	Not Reported
Al-Senani (4)	0.8±0.3°C/h	All patients	32.7±0.5°C
Arrich (5)	Not Reported	All patients	Not Reported
Pichon (43)	1.1±0.4°C/h	All patients	33±0.4°C in 91%
Kandzari (30)	3.3°C/h	All patients	33.55±0.35°C
Dixon (14)	Est: 3±0.85°C/h	All patients	Est: 33.75±0.5°C
DeGeorgia(12)	Not Reported	Not reached in 4/18	Not Reported
Georgiadis(20)	1.4±0.6°C/h	All patients	32.8±0.6°C
Guluma(22)	1.59°C/h	All patients	33.8±0.7°C
Hoedemaekers (58)	1.46°C/h: IVTM 0.72°C/h: SC	All 10 in IVTM 3/10 patients with SC	3.2% time out of range in IVTM (define as 0.2°C below or above target) 44 to 74% time out of range in SC
Merchant (59)	0.32°C/h: SC	Not reported	Target 32-34°C 63% < 32°C 28% < 31°C 13% < 30°C

IVTM: intravascular temperature management system; SC: surface cooling methods not included: Sunde(52)- didn't report these categories; Kliegel(32)- hypothermia induced with cold saline boluses

Central venous catheter-related complications

As clinicians consider the use of therapeutic hypothermia, they must carefully weight the risks as well as the potential benefits of the therapy. Medical complications associated with surface cooling techniques, such as shivering, thermal burns, limited patient access, difficulty during target temperature maintenance (overshoot), have been well described and shown in most cases not to be clinically significant. Endovascular cooling, however, involves the insertion of a central venous catheter into the femoral, internal jugular or subclavian vein, and the potential for central venous catheter-related complications. To be sure, this drawback is relative since most patients who would be considered candidates for therapeutic hypothermia require central venous access due to their medical condition. A recent retrospective study describing the use of femoral IVTM for 3-8 days in 10 severely brain-injured patients found subclinical (duplex ultrasound) evidence of catheter related thrombus in the 5 patients who had their catheters in place the longest.(50) None of these patients had a documented pulmonary embolism or any other clinically significant adverse consequence of the thrombus, however. (Of note, the manufacturer's recommendation is that the femoral cooling devices be left in place for no longer than 4 days-see Alsium IFU)

Some have suggested that subclinical thrombosis of the central venous system occurs in 20-30% of patients with indwelling catheters.(46) This is likely related to injury to the vessel endothelium from the catheter tip, thrombogenicity of the catheter itself, and length of time the catheter is left in the vein.(11) In addition, some venous insertion sites appear to be more prone to thrombosis than others.

Table 2: Reported Complications among Patients with Endovascular Cooling

Study	Disease	# of Subjects	Complications possibly related to catheters	Clinical or US Evidence of DVT
Feuchtl(17)	CA	19	Severe bradycardia req. pacemaker	NR
Sunde(52)	CA	29	Not Reported	NR
Flemming(18)	CA	31	Not Reported	NR
Al-Senani(4)	CA	13	Sepsis: possibly 2 cases;	No
Arrich(5)	CA	347: EC 114: SC	Not Reported	No
Pichon(43)	CA	40	Traumatic false aneurysm: 1 Infectious complications - 18 (no severe sepsis or shock): nosocomial bacteremia=5; pneumonia=5; bronchitis=6; UTI=2	No
Kandzari(30)	AMI	18	Not Reported	NR
Dixon(14)	AMI	20	Not Reported	NR
DeGeorgia(12)	Ischemic Stroke	18	Retroperitoneal hemorrhage - 1 Positive blood culture - 1 UTI - 1 Pneumonia - 2	3/18 (US)
Georgiadis(20)	Ischemic Stroke	6	Hiccups and pneumonia - all	No
Guluma(22)	Ischemic Stroke	10	No hematomas; No clinically evident vascular injury	No

Table 3: Reported DVT among Patients with Alsius Femoral Cooling Catheter

Study	Disease	# of Subjects	Possible DVT
Flemming(18)	CA	31	No
Sunde (52)	CA	29	No
Arrich (5)	CA	347	No
Feuchtl (17)	CA	19	No
Kliegel (60)	CA	9	No
Holzer (27)	CA	97	No
Kliegel (32)	CA	26	No
Al-Senani (4)	CA	13	No
Sunde (62)	CA	20	No
Pichon (40)	CA	40	No
Hoedemaekers C(58)	CA & TBI	10	No
Grimes (56)	CA	38	No
Thomas (61)	Stroke	10	No
Georgiadis (20)	Stroke	6	No
Keller (40)	SAH	90	No
Keller (31)	SAH	6	No
McDonagh (35)	Hypoxia	1	No
Megarbane (36)	Heat Stroke	1	No
Simosa HF (50)	TBI	10	Yes in 5 patients (Doppler US)
Lasater (38)	Accidental hypothermia	1	No

The reported incidence of catheter-related thrombosis depends on the detection method used, with clinical observation being the least sensitive, and direct ultrastructural examination of the catheter the most sensitive.(25) Duplex ultrasonography and phlebography are relatively sensitive techniques for the in vivo detection of clots, though the clinical significance of disease detected with these screening tools is unclear.

The most common insertion sites for central venous catheters are the subclavian, internal jugular, and femoral veins. In a study of 76 trauma patients who required unilateral femoral vein cannulation, and who had bilateral venous duplex sonography weekly for one month after cannulation, iliofemoral deep venous thrombosis was identified in 14% of the patients.(37) A study of 80 consecutive patients admitted to a mixed medical/surgical ICU found that the incidence of catheter related femoral vein thrombosis, as detected with phlebography, was 8.5%.(15) However, fibrin sleeves were observed in 15.7% of the cases. In another study of 124 mixed ICU patients who had femoral vein catheterization and were examined with duplex Doppler ultrasonography, 9.6% developed catheter related thromboses.(29) In that study, the incidence of thrombosis was not found to be related to duration of catheterization. Cumulatively, these studies suggest that the incidence of subclinical thrombosis with femoral central venous catheters is 8.5-26.2% in critical care patients with a wide variety of medical and surgical diseases.

To date there have been 21 peer-reviewed publications of the use of endovascular cooling catheters to treat patients with cardiac arrest, ischemic stroke, and several other diseases.(4-6, 12-14, 17, 18, 20, 22, 24, 30-32, 35, 36, 43, 45, 50-52) In patients who suffered cardiac arrest, acute myocardial infarction, ischemic stroke, 12 trials have been published. Six of those studies reported systematic evaluations for catheter related complications, and the incidence of such complications appears to be very low. (4, 12, 17, 20, 22, 43) (Table 2)

The DVT rate among all patients treated with the Icy catheter (Alsuis Corporation) also has been examined. Among 20 publications or presentations, only one of the studies reported duplex ultrasound evidence of thrombus, and then only in patients who had the catheter in place for more than 4 days. (Table 3) When all of the published clinical data is critically examined, the relative risk of Doppler ultrasound detected DVT is 0.6% (5/804).

Conclusions

IVTM is a new technology that, unlike surface cooling techniques, removes or transfers heat directly from the core thermal compartment. Compared with surface cooling methods, IVTM allows for a more efficient rate of cooling, tighter control at target temperature with the need for fewer nursing interventions, less shivering, and less thermal gradients between the peripheral and core compartment.

The benefits and risks for standard central venous catheters are well established. Lower extremity deep venous thrombosis is a known potential complication of a number of different acute medical and surgical diseases. Most critical care patients who would be candidates for therapeutic hypothermia will require central venous access to measure central venous pressure and deliver medications.

Multiple clinical studies have found that IVTM is associated with a less than 0.2% incidence of clinically significant thromboembolic disease. Among those with femoral catheters, the incidence of Alsuis endovascular cooling catheter-related thrombosis is very low. There is some evidence that the risk may be higher in trauma patients who have the catheter in place for more than 4 days. When considering methods to induce therapeutic hypothermia, or prevent fever, IVTM has many advantages compared to conventional cooling methods in terms of cooling effectiveness, and the ability to quickly reach and maintain the target temperature.

References

- Mild therapeutic hypothermia to improve the neurologic outcome after cardiac arrest. **N Engl J Med** 346:549-556, 2002.
- 2005 American Heart Association Guidelines for Cardiopulmonary Resuscitation and Emergency Cardiovascular Care. **Circulation** 112:IV1-203, 2005.
- Guidelines for the Management of Severe Traumatic Brain Injury, 3rd Edition. **J of Neurotrauma** 24, Supplement 1:S-26 - S-31, 2007.
- Al-Senani FM, Graffagnino C, Grotta JC, Saiki R, Wood D, Chung W, Palmer G, Collins KA: A prospective, multicenter pilot study to evaluate the feasibility and safety of using the CoolGard System and Icy catheter following cardiac arrest. **Resuscitation** 62:143-150, 2004.
- Arrich J: Clinical application of mild therapeutic hypothermia after cardiac arrest. **Crit Care Med** 35:1041-1047, 2007.
- Badjatia N, O'Donnell J, Baker JR, Huang D, Ayata C, Greer DM, Carter BS, Ogilvy CS, McDonald CT: Achieving normothermia in patients with febrile subarachnoid hemorrhage: feasibility and safety of a novel intravascular cooling catheter. **Neurocrit Care** 1:145-156, 2004.
- Bernard SA, Gray TW, Buist MD, Jones BM, Silvester W, Gutteridge G, Smith K: Treatment of comatose survivors of out-of-hospital cardiac arrest with induced hypothermia. **N Engl J Med** 346:557-563, 2002.
- Biswas AK, Bruce DA, Sklar FH, Bokovoy JL, Sommerauer JF: Treatment of acute traumatic brain injury in children with moderate hypothermia improves intracranial hypertension. **Crit Care Med** 30:2742-2751, 2002.
- Clifton GL, Allen S, Barradale P, Plenger P, Berry J, Koch S, Fletcher J, Hayes RL, Choi SC: A phase II study of moderate hypothermia in severe brain injury. **J Neurotrauma** 10:263-271; discussion 273, 1993.
- Clifton GL, Miller ER, Choi SC, Levin HS, McCauley S, Smith KR, Jr., Muizelaar JP, Wagner FC, Jr., Marion DW, Luerssen TG, Chesnut RM, Schwartz M: Lack of effect of induction of hypothermia after acute brain injury. **N Engl J Med** 344:556-563, 2001.
- Connors AF, Jr., Castele RJ, Farhat NZ, Tomaszewski JF, Jr.: Complications of right heart catheterization. A prospective autopsy study. **Chest** 88:567-572, 1985.
- De Georgia MA, Krieger DW, Abou-Chebl A, Devlin TG, Jaus M, Davis SM, Koroshetz WJ, Rordorf G, Warach S: Cooling for Acute Ischemic Brain Damage (COOL AID): a feasibility trial of endovascular cooling. **Neurology** 63:312-317, 2004.
- Diringner MN: Treatment of fever in the neurologic intensive care unit with a catheter-based heat exchange system. **Crit Care Med** 32:559-564, 2004.
- Dixon SR, Whitbourn RJ, Dae MW, Grube E, Sherman W, Schaefer GL, Jenkins JS, Baim DS, Gibbons RJ, Kuntz RE, Popma JJ, Nguyen TT, O'Neill WW: Induction of mild systemic hypothermia with endovascular cooling during primary percutaneous coronary intervention for acute myocardial infarction. **J Am Coll Cardiol** 40:1928-1934, 2002.
- Durbec O, Viviani X, Potie F, Viale R, Albanese J, Martin C: A prospective evaluation of the use of femoral venous catheters in critically ill adults. **Crit Care Med** 25:1986-1989, 1997.
- Eicher DJ, Wagner CL, Katikaneni LP, Hulsev TC, Bass WT, Kaufman DA, Horgan MJ, Languani S, Bhatia JJ, Givellichian LM, Sankaran K, Yager JJ: Moderate hypothermia in neonatal encephalopathy: efficacy outcomes. **Pediatr Neurol** 32:11-17, 2005.
- Feucht A, GB, Lawrenz T, Bartelsmeier M, Stellbrink C: Endovascular cooling improves neurological short-term outcome after prehospital cardiac arrest. **Intensivmed** 44:37-42, 2007.
- Flemming K, SG, Ziegs E, Diewok C, Gildemeister R, Wunderlich C, Strasser R: Comparison of external and intravascular cooling to induce hypothermia in patients after CPR. **GMS Ger Med Sci** 4, 2006.
- Gal R, Cundrie I, Zimova I, Smrcka M: Mild hypothermia therapy for patients with severe brain injury. **Clin Neurol Neurosurg** 104:318-321, 2002.
- Georgiadis D, Schwarz S, Kollmar R, Schwab S: Endovascular cooling for moderate hypothermia in patients with acute stroke: first results of a novel approach. **Stroke** 32:2550-2553, 2001.
- Gluckman PD, Wyatt JS, Azzopardi D, Ballard R, Edwards AD, Ferriero DM, Polin RA, Robertson CM, Thoresen M, Whitelaw A, Gunn AJ: Selective head cooling with mild systemic hypothermia after neonatal encephalopathy: multicentre randomised trial. **Lancet** 365:663-670, 2005.
- Guluma KZ, Hemmen TM, Olsen SE, Rapp KS, Lyden PD: A trial of therapeutic hypothermia via endovascular approach in awake patients with acute ischemic stroke: methodology. **Acad Emerg Med** 13:820-827, 2006.
- Higgins RD, Raju TN, Perlman J, Azzopardi DV, Blackmon LR, Clark RH, Edwards AD, Ferriero DM, Gluckman PD, Gunn AJ, Jacobs SE, Eicher DJ, Jobe AH, Laptook AR, LeBlanc MH, Palmer C, Shankaran S, Soll RF, Stark AR, Thoresen M, Wyatt J: Hypothermia and perinatal asphyxia: executive summary of the National Institute of Child Health and Human Development workshop. **J Pediatr** 148:170-175, 2006.
- Hinz J, Rosmus M, Popov A, Moerer O, Frensch I, Quintel M: Effectiveness of an intravascular cooling method compared with a conventional cooling technique in neurologic patients. **J Neurosurg Anesthesiol** 19:130-135, 2007.
- Hofbauer R, Moser D, Kaye AD, Dielacher C, Hornykewycz S, Handler S, Speiser W, Kapiotis S, Frass M: Thrombus formation on the balloon of heparin-bonded pulmonary artery catheters: an ultrastructural scanning electron microscope study. **Crit Care Med** 28:727-735, 2000.
- Holzer M, Behringer W, Janata A, Bayegan K, Schima H, Deckert Z, Losert U, Lagner AN, Sterz F: Extracorporeal venovenous cooling for induction of mild hypothermia in human-sized swine. **Crit Care Med** 33:1346-1350, 2005.
- Holzer M, Mullner M, Sterz F, Robak O, Kliegel A, Losert H, Sodeck G, Uray T, Zeiner A, Lagner AN: Efficacy and safety of endovascular cooling after cardiac arrest: cohort study and Bayesian approach. **Stroke** 37:1792-1797, 2006.
- Jiang J, Yu M, Zhu C: Effect of long-term mild hypothermia therapy in patients with severe traumatic brain injury: 1-year follow-up review of 87 cases. **J Neurosurg** 93:546-549, 2000.
- Joynt GM, Kew J, Gomersall CD, Leung VY, Liu EK: Deep venous thrombosis caused by femoral venous catheters in critically ill adult patients. **Chest** 117:178-183, 2000.
- Kandzari DE, Chu A, Brodie BR, Stuckey TA, Hermiller JB, Vetrovec GW, Hanman KL, Krucoff MW, Christenson RH, Gibbons RJ, Sigmon KN, Garg J, Hasselblad V, Collins K, Harrington RA, Berger PB, Chronos NA, Hochman JS, Calif RM: Feasibility of endovascular cooling as an adjunct to primary percutaneous coronary intervention (results of the LOWTEMP pilot study). **Am J Cardiol** 93:636-639, 2004.
- Keller E, Imhof HG, Gasser S, Terzic A, Yonekawa Y: Endovascular cooling with heat exchange catheters: a new method to induce and maintain hypothermia. **Intensive Care Med** 29:939-943, 2003.
- Kliegel A, Losert H, Sterz F, Kliegel M, Holzer M, Uray T, Domanovits H: Cold simple intravenous infusions preceding special endovascular cooling for faster induction of mild hypothermia after cardiac arrest—a feasibility study. **Resuscitation** 64:347-351, 2005.
- Kuboyama K, Safar P, Radovsky A, Tisherman SA, Stezoski SW, Alexander H: Delay in cooling negates the beneficial effect of mild resuscitative cerebral hypothermia after cardiac arrest in dogs: a prospective, randomized study. **Crit Care Med** 21:1348-1358, 1993.
- Marion DW, Penrod LE, Kelsey SF, Obrist WD, Kochanek PM, Palmer AM, Wisniewski SR, DeKosky ST: Treatment of traumatic brain injury with moderate hypothermia. **N Engl J Med** 336:540-546, 1997.
- McDonagh DL, Allen IN, Keifer JC, Warner DS: Induction of hypothermia after intraoperative hypoxic brain insult. **Anesth Analg** 103:180-181, table of contents, 2006.
- Megarbane B, Resiere D, Delahaye A, Baud FJ: Endovascular hypothermia for heat stroke: a case report. **Intensive Care Med** 30:170, 2004.
- Meredith JW, Young JS, O'Neil EA, Snow DC, Hansen KJ: Femoral catheters and deep venous thrombosis: a prospective evaluation with venous duplex sonography. **J Trauma** 35:187-190; discussion 190-181, 1993.
- Lasater et al. A novel method of intravascular temperature modulation to treat severe hypothermia. **Emerg. Med. J** 2006;23:56
- Metz C, Holzschuh M, Bein T, Woertgen C, Frey A, Frey J, Taeger K, Brawanski A: Moderate hypothermia in patients with severe head injury: cerebral and extracerebral effects. **J Neurosurg** 85:533-541, 1996.
- Keller E, et al. Experience with long-term hypothermia in 90 patients after aneurysmal subarachnoid hemorrhage. **Oral presentation at the 13th Annual International Symposium Intracranial Pressure (ICP) and Brain Monitoring, San Francisco, July, 2007**
- Nara J, Shiohagi T, Hara M, Saito I: Comparative effects of hypothermia, barbiturate, and osmotherapy for cerebral oxygen metabolism, intracranial pressure, and cerebral perfusion pressure in patients with severe head injury. **Acta Neurochir Suppl** 71:22-26, 1998.
- Nolan JP, Morley PT, Hoek TL, Hickey RW: Therapeutic hypothermia after cardiac arrest. An advisory statement by the Advancement Life support Task Force of the International Liaison committee on Resuscitation. **Resuscitation** 57:231-235, 2003.
- Pichon N, Amiel JB, Francois B, Dugard A, Etchecopar C, Vignon P: Efficacy and tolerance of mild induced hypothermia after out-of-hospital cardiac arrest using an endovascular cooling system. **Crit Care** 11:R71, 2007.
- Qiu WS, Liu WG, Shen H, Wang WM, Hang ZL, Zhang Y, Jiang SJ, Yang XF: Therapeutic effect of mild hypothermia on severe traumatic head injury. **Chin J Traumatol** 8:27-32, 2005.
- Schmutzhard E, Engelhardt K, Beer R, Brossner G, Pfäusler B, Spiss H, Unterberger I, Kampfl A: Safety and efficacy of a novel intravascular cooling device to control body temperature in neurologic intensive care patients: a prospective pilot study. **Crit Care Med** 30:2481-2488, 2002.
- Seigel EL, Jew AC, Delcore R, Iliopoulos JJ, Thomas JH: Thrombolytic therapy for catheter-related thrombosis. **Am J Surg** 166:716-718; discussion 718-719, 1993.
- Shankaran S, Laptook AR, Ehrenkranz RA, Tyson JE, McDonald SA, Donovan EF, Fanaroff AA, Poole WK, Wright LL, Higgins RD, Finer NN, Carlo WA, Duara S, Oh W, Cotten CM, Stevenson DK, Stoll BJ, Lemons JA, Guillet R, Jobe AH: Whole-body hypothermia for neonates with hypoxic-ischemic encephalopathy. **N Engl J Med** 353:1574-1584, 2005.
- Shiozaki T, Hayakata T, Taneda M, Nakajima Y, Hashiguchi N, Fujimi S, Nakamori Y, Tanaka H, Shimazu T, Sugimoto H: A multicenter prospective randomized controlled trial of the efficacy of mild hypothermia for severely head injured patients with low intracranial pressure. Mild Hypothermia Study Group in Japan. **J Neurosurg** 94:50-54, 2001.
- Shiozaki T, Sugimoto H, Taneda M, Yoshida H, Iwai A, Yoshioka T, Sugimoto T: Effect of mild hypothermia on uncontrollable intracranial hypertension after severe head injury. **J Neurosurg** 79:363-368, 1993.
- Simosa HF, Petersen DJ, Agarwal SK, Burke PA, Hirsch EF: Increased risk of deep venous thrombosis with endovascular cooling in patients with traumatic head injury. **Am Surg** 73:461-464, 2007.
- Steinberg GK, Ogilvy CS, Shuer LM, Connolly ES, Jr., Solomon RA, Lam A, Kassell NF, Baker CJ, Giannotta SL, Cockcroft KM, Bell-Stephens TE, Allgren RL: Comparison of endovascular and surface cooling during unruptured cerebral aneurysm repair. **Neurosurgery** 55:307-314; discussion 314-305, 2004.
- Sunde K, Pytte M, Jacobsen D, Mangschau A, Jensen LP, Smedsrud C, Draegni T, Steen PA: Implementation of a standardised treatment protocol for post resuscitation care after out-of-hospital cardiac arrest. **Resuscitation** 73:29-39, 2007.
- Takahashi I, Kitahara T, Endo M, Ohwada T: [Clinical analysis of hypothermia in children with severe head injury]. **No Shinkei Geka** 28:983-989, 2000.
- Tateishi A, Soejima Y, Taira Y, Nakashima K, Fujisawa H, Tsuchida E, Maekawa T, Ito H: Feasibility of the titration method of mild hypothermia in severely head-injured patients with intracranial hypertension. **Neurosurgery** 42:1065-1069; discussion 1069-1070, 1998.
- Tokutomi T, Morimoto K, Miyagi T, Yamaguchi S, Ishikawa K, Shigemori M: Optimal temperature for the management of severe traumatic brain injury: effect of hypothermia on intracranial pressure, systemic and intracranial hemodynamics, and metabolism. **Neurosurgery** 52:102-111; discussion 111-102, 2003.
- Grimes Anderson, Horuchi, Concha, Fleegler, Hurwitz et al. Resuscitative Hypothermia After Cardiac Arrest: Performance in a Community Hospital. **CHEST / 128 / 4 / October, 2005 Supplement 167S.**
- Zhi D, Zhang S, Lin X: Study on therapeutic mechanism and clinical effect of mild hypothermia in patients with severe head injury. **Surg Neurol** 59:381-385, 2003.
- Hoedemaekers CC, EzzahatMM, et al: Comparison of different cooling methods to induce and maintain normo- and hypothermia in ICU patients: a prospective intervention study. **Critical Care**. 2007,11:R91
- Merchant RM, Abella BS, et al. Therapeutic hypothermia after cardiac arrest: Unintentional overcooling is common using ice packs and conventional cooling blankets. **Crit Care Med** 2006 Vol. 34, No. 12
- Kliegel et al : Cold infusion alone are effective for induction of therapeutic hypothermia but do not keep patients cool after cardiac arrest. **Resuscitation** (2007) 73, 46 - 53
- Thomas et al. Safety and Therapeutic Benefit of Hemicraniectomy Combined with Mild Hypothermia in Comparison with Hemicraniectomy Alone in Patients with Malignant Ischemic Stroke; **Cerebrovascular Diseases** 2006;21:79-85
- Sunde K, et al. Therapeutic hypothermia with Endovascular cooling. **Intensive Care Med**. 2004;12:1



15770 Laguna Canyon Road, Suite 150
Irvine, California 92618-3111

USA
Tel: +1-949-453-0150
Fax: +1-949-453-0702

www.alsius.com

MARION0907