

Masimo SET

Green indicator light  
Red triangle warning icon  
Yellow triangle warning icon  
OK button

**A  
E  
V**

HR 120  
stby 40  
SpO2  
stby 94  
FiO2 21  
O2 USE 0.0  
PIP 10 35  
PEEP 5 off  
Vt<sub>(del)</sub> 1567 off off  
BPM 14 40 4  
Mode NPPV  
CPAP

CPAP+NPPV+PS(2)

Vmin 32.5 MAP 6.3 TI 1.14

100%

35  
10  
5

↓ MENU ↓

Menu button

Cancel button (X)

Confirm button (✓)

Manual Breath button

**ZOLL**

# The Most Prehospital-Ready Ventilator™ Ready for Your Crews

## Smart Help™

Smart Help guides the user with prompts about the nature of the alarm and how best to manage it.

## Apnea Backup

ZOLL ventilators deliver pressure-targeted breaths if a spontaneous breath is not detected when operating in a therapeutic mode.

## Automatic Leak Compensation

ZOLL ventilators automatically detect an ill-fitting mask or leaking airway and provide extra flow when operating in a therapeutic mode.



## Touch

Touch the button that corresponds to the parameter you want to change.

## Turn

Turn the dial to adjust that parameter to the desired value.

## Confirm

Confirm this new value by pressing the Confirm/Select button.

# ZOLL®

To learn more about the AEV ventilator, visit us at  
[www.zoll.com](http://www.zoll.com) or call us at 1-800-804-4356.